EN

VACUUM EXTRACTOR MODIFIED BY BIRD

VACUUM EXTRACTOR MODIFIED BY BIRD

This obstetrical instrument is utilized to assist a clinician in the delivery of an infant during childbirth.

Important note

Proper surgical procedures and techniques are the responsibility of the medical practitioner. Each physician must evaluate the appropriateness of the procedure based on his or her own medical training and experience.

General information

The Medela stainless steel vacuum extractors come in two different models. The original model by Malmström and the modified model by Bird. Each model consists of three stainless steel suction cups in the sizes 40, 50 and 60 mm diam. The cups of each model are identical, but the mode of traction and tube connectors vary

Pump requirements

- This obstetrical instrument has been developed for use in conjunction with an external vacuum source capable of precise vacuum regulation to ensure proper vacuum levels (up to -80 kPa/-600mmHg).
- The vacuum source should be designed for medical purposes and may be either electric or manual.
- We recommend using the Medela Basic suction pump.
- Note: Wall suction is not to be used.**

Tubing requirements

The tubing connecting the instrument to the vacuum source is vital for safe and proper function of the instrument. To ensure this only Medela original tubing should be used.

Cleaning/Disinfection/before first use

- After use, soak the suction cups and tubings at room temperature for 30 minutes in disinfect-tion/cleaning solution A).
- After the reaction time, rinse the residue of the solution from the parts with water. Rinse in clean tap water.
- Remove locking ring, tubing and bottom plate (with a finger, with forceps or by pushing a small wooden rod through the tube-connector-piece of the cup).
- Clean all metal parts, locking ring and bottom plate with a brush and cleaning agent B).
- Replace the bottom-plate with a new one when satisfactory cleanliness cannot be achieved.
- Rinse the tubing for at least 2 minutes with the cleaning solution B).
- The parts should be free from residue and foreign particles.
- Rinse the parts and the tubing with water (preferably fully desalinated) for 2 minutes.
- Then dry the parts with a fiberless cloth.
- Place the polypropylene bottom plate in position inside the cup.
- A) Fluid aldehyde-free disinfection solution with cleaning effect for use as soaking bath, bactericide, fungicide, limited virucide, with good material compatibility for stainless steel, non-ferrous material and plastics including silicone, slightly alkalisch.
- B) Fluid, pH-neutral, enzymatic cleaning solution for the treatment of instruments by machine or manually with very good material compatibility for stainless steel, non-ferrous material and plastics including silicone.

For the validation of the procedure a 1% neodisher® Med2Zm-solution was used.

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• The regulations of the manufacturer of the cleaning agent, resp. use and dangers must be kept.

• The use of other cleaning and disinfecting agents should be avoided with the hospital internal hygiene administrator.

• **Sterilisation (also before first use)**

- Assemble and pack the suction cups before sterilisation.
- Place the polypropylene bottom plate in position inside the cup.
- Connect the long silicone rubber suction tubing to the tube-connector piece of the cup and secure it by means of the locking ring.
- Sterilisation
 - In the autoclave with saturated steam at 134°C, with triple pre-evacuation, sterilisation time 6 minutes.

Expected product life

- For the metal cup, chain and handle:
 - 5 years.
- Tubing:
 - 60 cleaning/sterilisation cycles.

Safety rules

- If damages - especially to the chains and tubings (tears) are determined, the parts may no longer be used.
- No modification of this equipment is allowed.
- Do not use it:
 - If the device is damaged.
 - If the device shows clear safety defects.
- Wall suction is not to be used.
- Vacuum extraction must only be performed or supervised by trained and experienced clinicians.
- The neonatal care staff must be aware of the indications, contraindications and precautions for assisted delivery and must watch for the signs of complications of vacuum assisted delivery.

- Operative vaginal delivery should not be performed when the probability of success is very low.
- Ensure against intrusion of extraneous tissues and ensure proper placement of cup.
- Attention: Never apply the cup to any portion of the infant's face as this can cause fetal trauma.
- The neonatal care staff must be aware of the indications, contraindications and precautions to avoid malfunctions, maternal or fetal trauma and maternal discomfort.
- The head, not just the scalp, must descend with each pull.
- Pull only during a contraction.
- Always apply traction at an angle allowing the stem of the cup to remain in-line or on-axis.
- Never apply rocking movements or torque to the device; only steady traction should be used.
- Discontinue vacuum extraction if there is evidence of significant scalp trauma or damage.
- The head must be delivered, or almost completely delivered, with no more than three pulls.
- The head must be completely delivered within 15–20 minutes of first applying the cup.
- Never exceed 10 minutes of cumulative traction time.
- Never exceed maximum recommended vacuum -80 kPa (-600 mmHg) or "pop offs" (2) - scalp trauma can occur.
- Discontinue vacuum extraction if the vacuum cup becomes disengaged or pops off twice.
- Do not pull three times, if the second pull fails to produce good descent.
- Parents should be instructed to check the baby's health good.
- The infant has to be watched for several days for possible problems.

Disposal

Dispose of according to local regulations.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

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